

# CES Device Use Cases

## **Team 23:**

- *David Chen*
- *Zachary Hamel*
- *Minh Thang Cao*

## **UC 1: Turn Device On**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The Device is Off

Success Guarantee:

- The Device turns On

Main Success Scenario:

- 1) The Patient presses the Power Button
- 2) The Device Turns On

Extensions:

2a. Battery runs out: The Device Does Not Turn On

2b. The device disabled itself permanently due to exceeding 700 amperes: The Device Does Not Turn On

## **UC 2: Turn Device Off**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The Device is already on

Success Guarantee:

- The Device turns off

Main Success Scenario:

- 1) User Presses the power Button
- 2) The Device Turns off

## **UC 3: Administer Treatment**

Primary actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- User already attached the earclips on their earlobes and chose the frequency based on medical needs

Success Guarantee:

- The treatment is completed

Main Success Scenario:

- 1) User turns on the device (include UC 1)
- 2) The Device's Battery begins to drain (include UC 7)
- 3) User sets the Timer (20, 40, 60 minutes) (include UC 5)
- 4) User sets the Treatment Frequency (include UC 4)
- 5) User sets the Treatment Waveform (include UC 10)
- 6) The User presses the Start Button to start the treatment
- 7) The Device circulates the treatment current through the user's head
- 8) User changes the treatment current multiple time based on how they feel during the process (include UC 6)
- 9) The Device turns off automatically at the end of timed cycle

Extensions:

- 1a) The device is out of power and does not turn on
- 1b) When the device is on, the device loses power and turns off before the treatment is complete
- 3a) User doesn't have earclips on (electrodes are not in contact with skin): Device pauses the process until the electrodes are in contact with the skin
- 3b) Instead of choosing the treatment time, the User decides to play a recorded treatment (refer to UC 8)
- 4a) The treatment exceeds 700 millimeters and automatically disables itself for the users' safety
- 5a) The User does not press the start button for 30 minutes and the device shuts off.
- 6a) The User clicks a Timer option, resetting the Timer of the device to the selected option
- 6b) The earclips fall off the user. If skin contact is lost during treatment for less than 5 seconds treatment resumes, otherwise, treatment stops.
- 9a) Before the device turns off, the user decides to record the treatment (refer to UC 9)

## **UC 4: Change Frequency**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The device is on

Success Guarantee:

- Frequency is changed for the current treatment

Main Success Scenario

- 1) Patient can choose between 0.5 Hz, 77Hz and 100 Hz by clicking frequency button
- 2) Frequency of the current treatment is changed to the selected value

Extensions:

- 1a) User choose the selected frequency: nothing will be changed

## **UC 5: Choose Countdown Time**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The treatment has not started yet

Success Guarantee:

- Countdown timer is set for the current treatment

Main Success Scenario:

- 1) Patient can choose between 20 minute, 40 minute, and 60 minute treatment times
- 2) Patient confirms by pressing the Start Button
- 3) The Device's treatment will last for the specified time

## **UC 6: Change Treatment Current**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The treatment has started

Success Guarantee:

- Treatment electrode current is changed successfully as input

Main Success Scenario:

- 1) If a patient wants to increase treatment level, the patient presses the up button. If the patient wants to decrease treatment current, the patient presses the down button.
- 2) The device current either increases by 50 amperes or decreases by 100 amperes
- 3) Patient can continue increasing or decreasing until the current is at 500 or 0 amperes.

Extensions:

3a. The treatment exceeds 700 millimeters and automatically disables itself for the users' safety

## **UC 7: Change Battery Level**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The device is on

Success Guarantee:

- Treatment electrode current is changed successfully as input

Main Success Scenario:

- 1) The device is turned on
- 2) The device drains battery based on the time used

Extension:

- 2a. User increase treatment current: the device drains battery faster based on the treatment current set
- 2b. If the device reaches battery level 5% the device will issue a warning to the user
- 2c. If the device reaches battery level 2% the device will automatically shut off

## **UC 8: Record a Therapy**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- User is just about to start a treatment

Success Guarantee:

- The treatment is recorded

Main Success Scenario

- 1) User check the Record check box
- 2) The user then administers a treatment
- 3) The Treatment is then recorded in the device

## **UC 9: Replay a Recorded Therapy**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- User has recorded at least one therapy
- User has earclips on

Success Guarantee:

- The therapy process is done exactly the same as recorded

Main Success Scenario:

- 1) User chooses one of the treatment in the history list
- 2) Device execute a treatment with the same settings

## **UC 10: Change Waveform**

Primary Actor: The Patient

Level: User Goal

Scope: The Device

Precondition:

- The device is on

Success Guarantee:

- Waveform is changed for the current treatment

Main Success Scenario:

- 1) Patient can choose between Alpha, Beta or Gamma by choose one of the radio buttons
- 2) The current treatment's waveform is changed to to selected one

Extensions:

- 1a) User choose the selected waveform: nothing will be changed